




Contralateral Foraminal Area Increases Significantly After MIS Transforaminal Lumbar Interbody Fusion Using Biplanar Expandable Cage

Yee Gen Lim, David Shaoen Sim, Arellano Pasion Dizon,
Chuen Jye Yeoh, Lei Jiang, Zhixing Marcus Lim

Department of Orthopaedic Surgery, Singapore General Hospital

Introduction

- Biplanar expandable (BE) cages recently designed for use in transforaminal lumbar interbody fusion (TLIF)
- 2-piece device consisting of a polyetheretherketone (PEEK) shell that expands bidirectionally with the insertion of a titanium shim
- Open-architecture design which allows backfilling of bone graft and conforms to patient's endplate configuration
- Limitations of current designs:
 - Static cages – small surgical window
 - Uniplanar expandable cage – small contact area risks endplate violation and cage subsidence
- Advantage of BE cages – expansion in width and height
 - Improved restoration of disc height and segmental lordosis
 - Decreased theoretical risk of subsidence due to greater implant-endplate contact surface area



Introduction

- MIS-TLIF performed through unilateral laminectomy and facetectomy on more symptomatic side
- Management of contralateral neuroforaminal foraminal stenosis depends on indirect decompression resulting from restoration of disc height and listhesis correction.
- BE cages have been shown to improve early to mid-term clinical outcomes and radiographic parameters
- However radiographic parameters were only measured with plain radiographs and unable to measure foraminal parameters




1. Sim, D. S., Kasivishvanaath, A., Jiang, L., Cheong Soh, R. C., & Ling, Z. M. (2023). Biplanar expandable cages for transforaminal lumbar interbody fusion are safe and achieve good 1-year clinical and radiological outcomes in an Asian population. *International Journal of Spine Surgery*, 17(4), 520–525. <https://doi.org/10.14444/8472>
2. Coric D, Roybal RR, Grubb M, et al. Bidirectional expandable technology for transforaminal or posterior lumbar interbody fusion: a retrospective analysis of safety and performance. *Int J Spine Surg*. 2020;14(s3):S22–S30. doi:10.14444/7123
3. Tan LA, Rivera J, Tan XA, Le VP, Khoo LT, Berven SH. Clinical and radiographic outcomes after minimally invasive transforaminal lumbar interbody fusion-early experience using a biplanar expandable cage for lumbar spondylolisthesis. *Int J Spine Surg*.2020;14(s3):S39–S44. doi:10.14444/7125


Objective/Aim




- Investigate change in radiographic parameters including neuroforaminal parameters after MIS-TLIF with biplanar expandable cages

Methods

- Retrospective study on consecutive patients who underwent MIS-TLIF utilizing biplanar expandable cages in 2023 from a single institution
- Intraoperative 3D imaging performed before and after MIS-TLIF were analyzed to calculate change in radiographic parameters
- Inclusion criteria:
 - Degenerative disc disease, spondylolisthesis or spinal stenosis
 - 1-level, 2-level and 3-level lumbar fusions
 - Minimally invasive (MIS) technique
- Exclusion criteria: Previous spinal fusion, infection or malignancy
- Radiographic parameters assessed:
 - Ipsilateral and contralateral foraminal area
 - Posterior and anterior disc height, listhesis, segmental lordosis,







Demographics/Patient Characteristics

Table 1. Patient Characteristics	
No. of patients	12
Age (years)	70.8 ± 5.8
BMI	26.5 ± 3.3
Gender	3 male and 9 female
Levels	5 1-level TLIF and 7 2-level TLIF
	19 spinal levels fused

Results

- Significant improvement in posterior disc height ($+2.7 \pm 1.2 \text{mm}$; $p < 0.001$), anterior disc height ($+2.6 \pm 3.1 \text{mm}$; $p < 0.001$), listhesis, segmental lordosis ($+3.2 \pm 2.6^\circ$, $p < 0.001$)
- Significant improvement in ipsilateral ($41.3 \pm 17.5 \text{cm}^2$; $p < 0.001$) and contralateral foraminal area ($29.1 \pm 16.7 \text{cm}^2$; $p < 0.001$)
- Ipsilateral and contralateral foraminal area increased by 41% and 31% respectively

Table 2. Radiological Outcomes				
BE cage (n=12)	Pre-op	Post-op*	Difference	p-value
Radiological Outcomes				
Segmental Lordosis($^\circ$)	6.7 ± 4.0	10.0 ± 3.4	3.2 ± 2.6	$< .001$
Anterior Disc Height(mm)	8.4 ± 2.1	11.0 ± 3.1	2.6 ± 3.1	$< .001$
Posterior Disc Height(mm)	4.3 ± 1.5	7.0 ± 1.4	2.7 ± 1.2	$< .001$
Foraminal Area(cm^2)				
- Ipsilateral	1.1 ± 0.3	1.5 ± 0.2	0.4 ± 0.2	$< .001$
- Contralateral	1.0 ± 0.2	1.3 ± 0.2	0.3 ± 0.2	$< .001$

Conclusion

- Utilizing biplanar expandable cage in MIS-TLIF can lead to significant improvement in radiographic parameters, including increase in contralateral foraminal area